



Tulsi / Basil

Ocimum basilicum

ADAPTOGEN

Botanical / Latin name	Ocimum basilicum
Category	Adaptogen
Standardisation	Marker-standardised extract - Typical 1%-5% active (HPLC). Available up to 10%+ for select markers.
Marker compound	Ursolic acid (method: HPLC)
Available grades	2% ursolic acid; 10:1 ratio
Available forms	Dry, Liquid
Solubility & grades	Dry extract: water-soluble & hydro-ethanolic profiles. Spray-dried and granulated grades available.
Applications	Stress & cortisol support; sports recovery; sleep & adaptogenic stacks; functional beverages; premium capsule lines.
Appearance	Fine powder (typical); varies by form
MOQ	On request
Shelf life	24 months (typical)
Packaging	HDPE drums, food-grade liners
Lead time	2-4 weeks
Certifications	GMP Certified, ISO 9001:2015, NABL tested, FSSAI compliant
Typical potency	Typical 2% (range 1-2.5%) confidence: low

Custom standardisation, ratios, particle size and packaging available to your specification.

Request a quote or a sample with batch COA - we reply within 24 hours.